

III. Physikalisches
Institut A

RWTHAACHEN
UNIVERSITY


Experimental Techniques in Particle Physics (WS 2020/2021)

Exercises: the CMS muon DT system

08.12.2020


Prof. Alexander Schmidt

Intro




Intro

schematic of a DT layer
consisting of many cells



Intro



one DT cell:



Important facts:

- wire diameter 50 μm
- 85% / 15% Argon-CO₂ gas mixture
- designed in a rectangular way with almost homogeneous field !!
- the drift velocity is therefore almost constant $v=54 \mu\text{m}/\text{ns}$
- the maximum drift time is therefore 390 ns (15 bunch crossings !!!)
- muons go through almost straight upright
- reconstruction of the **position** of the charge deposit possible from the drift time t_{drift} (this leads to a reconstructed hit position \rightarrow "rechit")
- left-right ambiguity is resolved through staggered arrangement

Intro:



- the drift time (for electrons) is read out by the CMS data acquisition system
- it needs to be calibrated with other time effects:
 - propagation time of the signal along the anode wire
 - cable length
 - read-out electronics processing time
 - trigger latency

Exercise:

- A file is provided in jupyter in the directory ex6-CMS_DT
- the file contains 1000 drift times (in ns) which are already pre-calibrated
- the file also contains the true rechit position (in cm), which is an extrapolation from a track using other layers (or true position based on simulation)
- from the given information: calculate the **resolution** of the CMS muon rechits
- compare your result with the official documentation available from the CMS collaboration