



III. Physikalisches
Institut A

RWTHAACHEN
UNIVERSITY

Experimental Techniques in Particle Physics (WS 2020/2021)

Exercises

Prof. Alexander Schmidt

01.12.2020

Exercise

draw the time-dependent signal (current and/or voltage) created in a gas detector for the following cases:

- **parallel plate drift chamber**
 - distance of plates 1cm
 - size of plates 1m²
 - e⁻-I⁺ pair deposited in the middle between the two plates
 - Argon gas at atmospheric pressure
 - Voltage between plates 2kV
- **cylindrical proportional chamber**
 - tube radius 1cm
 - wire radius 10μm
 - 100 e⁻-I⁺ pairs deposited in the middle between wire and tube
 - same gas and Voltage as above
 - hint: electron drift is negligible compared to avalanche and resulting ion drift, the gas Gain factor G=10000
- **for next week:**
 - how does the signal of a continuous track look like
 - the track is created by a 1 GeV muon (calculate the ionisation and the drift of the ions)